

# Tilty Vario



*Tilttable wheelchair for  
adolescent and adults  
Active, tilttable and growable*



# Tilty Vario

Active, tiltable, growable

German aid index N° 18.99.02.1013



seat width  
**28-44 cm**  
(+2 cm)



back height  
**30-50 cm**  
(+5 cm)



seat depth  
**30-44 cm**  
(+4 cm)



weight capacity  
**100 kg**



wheel sizes  
**22" and 24"**




tilting angle  
from **0° to +25°**



## Equipment examples

All equipment options are on the order form of the Tilty Vario



holder to attach  
abduction wedge




truss pads



arm pads adjustable in  
height and angle



undivided foot rest



head rests with chest-  
shoulder-belt

And  
much  
more ...

## Overview

- Tilttable wheelchair for adolescents and adults
- Back adjustable in angle (80° – 120°) and completely stowable
- Highly stable rigid frame with 10 cm abduction; also suitable for Pörnbacher therapy methods
- Smooth handling with minimized effort
- Suitable as a seat for transport in a vehicle according to ISO 7176-19